

Comfort-3

Eco

THERAPY CHAIR

HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY

- The Comfort-3 ECO Therapy Chairs are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably and therefore help to relieve tension during treatment
- Several optional accessories for optimal choice
- Slim design, less space needed



PROFESSIONAL SEATING SOLUTIONS BY

Digiterm[®]
medical equipment

* Chair in the photo involves some optional elements

Comfort-3 Eco

ART.-NR. 723C3E00

Optional accessories	Art.-Nr.
EPN kit incl. multicontact pin with internal connecting cable HIGHLY RECOMMENDED FOR DIALYSIS	713E004
Trendelenburg foot-switches on both sides of underframe	713E005
Rechargeable battery back-up operation in case power failure (24V DC)	713E006
Reading lamp with LED light	713E008
USB charging unit for smart devices (5V, 1A)	713E009
Connecting cable (3 m) for potential equalization (EPN)	713E012
Comfort-T2 Multifunctional table	713M009
Flexible height adjustable armrest for universal medical applications (one pair)	713M010
Integrated gas spring in the arm support for smoother and safer operation (one pair)	713M011
Push bar	713M012
Universal holder for tablets and smartphones	713M013
Handle on both sides of the seat	713M014
Motorized height adjustment	713M017
Central operating brake with 4x diam. 7.5 cm twin wheel castors	713M020
4x diameter 10 cm swivel castors with total lock	713M021
Paper roll holder	713M033
Infusion bottle holder (IV pole) with height adjustment	713M038
Individual logo on the upholstery	713U011
Headrest pillow square	713U012
DuoShape upholstery in individual colour combination (40 diff. colours) TUNDRA	713U026
Transparent protecting plastic cover for short legrest (3 segments upholstery)	713U029
3 segments upholstery in individual colour combination TUNDRA	713U040

General technical data			
Fixed seat height	55 cm	Size of castors	4 x Ø 7.5 cm twin wheel castors with total lock
Optional height adjustment range	55–70 cm	Weight of Chair	60 kg
Backrest adjusting range	80° (-14° under horizontal in Trendelenburg position)	Maximal load	200 kg weight capacity
Seat tilt adjustment	14°	Standard Power Supply	100-240 V AC, 50/60 Hz
Legrest adjustment	90°	Technical data of motors	24V DC, IPX6
Dimensions in seat position	L 99 x W 92.5 x H 139 cm	Frame colour	RAL 9001
Size of upholstery	L 195 x W 60 x H 9.5 cm	Applied directive and standards	93/42/EEC and harmonized standards (EN60601-1; EN60601-2-52; etc.)
Armrest size	L 57 x W 15 x H 6.5 cm		



Digiterm Ltd.
Ipar út 7, H-3200 Gyöngyös
Tel.: +36 37 311 990
Fax: +36 37 311 893
E-mail: sales@digiterm.hu
www.digiterm.hu



Tested Medical Device Safety Regular Production Surveillance



www.tuv.com
ID 1111237221

ISO 9001
ISO 13485
ISO 14001