

***Baxter***

## Welch Allyn Stethoscopes

Hear the Difference.



# Welch Allyn Stethoscopes

Hear the difference from a leader in diagnostic care.

From our **Harvey DLX** cardiology stethoscope with trumpet brass construction designed to enhance sound, to our well-known professional line, every **Welch Allyn** product has been designed with you — the clinician — in mind.



## Cardiology Grade Stethoscopes

Developed during a 30-year association with Dr. W. Proctor Harvey, MD, a recognized “Master” by the American College of Cardiology.

### Welch Allyn Harvey DLX Stethoscope

The trumpet brass construction is designed to help you hear more of the important cardiac, pulmonary, and vascular sounds without relying on “variable pressure” techniques.

### Welch Allyn Harvey Elite Stethoscope

Engineered with the exceptional comfort and auscultation capabilities physicians have come to rely on. Its stainless steel head has the ability to identify the lowest frequencies detectable by the human ear.



## Welch Allyn Professional Grade Stethoscopes

Durability and comfort for the practicing professional — they’re lightweight and can be custom fitted via adjustable binaurals and eartips — designed for listening to cardiac, pulmonary and vascular sounds.

### Professional

- Double-head design
- Flat diaphragm
- Adult and Pediatric versions available in up to six colors
- Stainless steel chestpiece and binaurals
- Rotatable binaurals and interchangeable comfort sealing eartips
- Not made with natural rubber latex
- 5-year warranty

### Lightweight

- Chrome-plated chestpiece for durability
- Rotatable binaurals
- Double-head design
- Flat diaphragm
- Adult version only available in four colors
- 2-year warranty





## Welch Allyn Stethoscope Specifications

	Cardiology Grade		Professional Grade		
Product name	Harvey DLX (28")	Harvey Elite Double Head (22") (28")	Professional Adult	Professional Pediatric	Lightweight
Product numbers	5079-321 Triple, Black 5079-322 Triple, Burgundy 5079-323 Triple, Navy 5079-324 Triple, Green 5079-325 Double, Black 5079-326 Double, Burgundy 5079-327 Double, Navy 5079-328 Double, Green	5079-122 Elite, Black (22") 5079-125 Elite, Black 5079-270 Elite, Burgundy 5079-271 Elite, Navy 5079-284 Elite, Green 5079-125P Elite, Ped Blk	5079-135 Black 5079-137 Navy 5079-138 Plum 5079-139 Burgundy 5079-285 Forest 5079-289 Blue	5079-145 Black 5079-147 Navy 5079-148 Plum 5079-149 Burgundy 5079-287 Forest	5079-73 Velvet Black 5079-74 Poppy Red 5079-75 Dawn Gray 5079-76 Misty Blue
Product Type	Double/Triple Head	Double Head	Double Head	Double Head	Double Head
Head Diameter Inches	1.9 in.	1.9 in.	1.75 in.	1.2 in.	1.9 in.
Overall Length Inches	28 in.	28 in.	28 in.	28 in.	32 in.
Tubing Length Inches	17 in.	17 in.	21.5 in.	21.5 in.	24 in.
Chestpiece Material	Brass	Stainless Steel	Stainless Steel	Stainless Steel	Chrome Plated Zinc
Diaphragm Material	Polycarbonate	Polycarbonate	Clear Polycarbonate	Iten Resiten G-10/FR4T Glass Epoxy Sheet	PVC
Diaphragm Style	Flat or Corrugated	Flat or Corrugated	Flat	Flat	Flat
Adjustable Binaural	Yes	Yes	Yes	Yes	Yes
Eartips	Soft & Hard	Soft & Hard	Soft & Hard	Soft & Hard	Hard
Tubing Type	Dual Lumen	Dual Lumen	Single Lumen	Single Lumen	Single Lumen
Anti-Chill Ring	Yes (Brass)	Yes	Yes	Yes	Yes
Ships Assembled	Yes	Yes	Yes	Yes	No
Packaging	Printed Box	Printed Box	Printed Box	Printed Box	Plain White Box 1-piece
Description	Designed to meet the demanding needs of cardiologists, this stethoscope has fully adjustable binaurals and triple leaf spring for comfort. The head is made of trumpet brass because of its ability to transmit sounds accurately. The dual-lumen tubing is not made with natural rubber latex.		Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs. Stainless-steel chestpiece intensifies body sounds transmitted through tissues. Comfort sealing eartips depress and conform to the shape of ear canals for soft yet airtight fit. Not made with natural rubber latex.		Quality design for blood pressure measurement and general auscultation. Lightweight and flexible single-channel tubing provides comfortable use over extended periods. Adjustable chrome-plated binaurals.

**Rx only.** For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

### Baxter.com

Baxter International Inc.

Baxter, Harvey DLX, Harvey Elite, Hillrom and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC137-240002 (v1.0) 05/2025