

PRODUCT CATALOGUE

EVERYTHING YOU NEED AT THE POINT OF CARE



WelchAllyn®

SOLUTIONS FOR PRIMARY CARE AND CLINICS

Your primary care partner, helping you:

- Perform more comprehensive examinations—faster
- Better diagnose and manage chronic conditions
- Improve financial outcomes



At Hill-Rom, we're passionate about helping caregivers address the most pressing challenges they face every day.

Our portfolio of solutions includes technologies that help your staff reduce preventable errors, programmes that reduce falls and improve patient mobility and best-in-class products that provide the accurate data clinicians need to anticipate care.

It's all designed to help make your staff more efficient, reduce risk and keep your patients safe.

One Hill-Rom. Your premier partner to enhance outcomes for patients and their caregivers.

WelchAllyn®

Hill-Rom®

Pocket LED Oscopes, Ophthalmoscopes and Sets

Maximised performance and durability in a compact set that travels easily with you from one patient to the next.

- > Examine your patients in style with high-performance diagnostic tools
- > Available in two versions: Pocket LED and Pocket Plus LED with enhanced features
- > Protect and customise them with handle bumpers and otoscope window bumpers to mix and match your look
- > Lightest-weight Welch Allyn Ophthalmoscope and Otoscope for portability and ease of use
- > The most durable, pocket-style otoscope on the market today
- > Featuring Welch Allyn SureColour LED bulb technology

Choose the set that's right for you

	Pocket LED	Pocket Plus LED
Light output	3 lumens	3 lumens (ophth) 5.5 lumens (oto)
Ophth neutral density filter	No	Yes
Accessories	2 sets	4 sets
Case	Soft	Hard
Warranty	3 years	5 years



- 92871-BLK** Pocket LED Set - Black with Soft Case
- 92871-WHT** Pocket LED Set - White with Soft Case
- 92871-PUR** Pocket LED Set - Plum with Soft Case
- 92871-BLU** Pocket LED Set - Blueberry with Soft Case
- 92880-BLK** Pocket Plus LED Set - Black with Hard Case
- 92880-WHT** Pocket Plus LED Set - White with Hard Case
- 92880-PUR** Pocket Plus LED Set - Plum with Hard Case
- 92880-BLU** Pocket Plus LED Set - Blueberry with Hard Case

Individual otoscopes and ophthalmoscopes also available, please call customer care for more information.


92871-BLU



92880-PUR




OPTIMISED LED



WELCH ALLYN
SURECOLOUR™ LED


30x LONGER LIFE


- > Tested to 30x halogen bulbs
- > Guaranteed for 7 years



ENERGY EFFICIENCY

- > 3x the on-time vs. halogen with Welch Allyn 3.5V rechargeable handles
- > Lower operating temperature intended for greater patient/clinician comfort





For more information about our products and solutions please visit www.welchallyn.co.uk or contact customer care on 0207 365 6780

WelchAllyn®